


# Waters Xevo G2-xs Quadrupole Time-of-Flight (Q-TOF MS)

## University of Alabama Mass Spectrometry Facility



**ESI Source**

**Capillary:** 2-4 kV

**Sample Cone:** 25-100 V

**Source Offset:** 25-100 V

**Source temperature:** for the source block, 60-120°C

**Desolvation gas :** 300-1000 L/hr

**Desolvation gas temperature:** 150-650°C (typically 150-300°C), depends on solvent


**Cone gas flow rate:** 0-150 L/hr

**RF setting**


StepWave: 300 V

Ion guide: 350 V

**StepWave 2 offset:** 25 V



- Analytes enter the large diameter ion guide
- Ions enter the upper ion guide (by using differential voltage); neutrals and solvents are directed to the pump inlet




**Quadrupole**

Arbitrary value calculated from the RF/DC ratio

LM Resolution	4.7 (4k quad) 12 (32k quad)	Sets the width of the quad transmission window in Quad and TOF MS/MS modes.
HM Resolution	15	Sets the width of the quad transmission window in Quad and TOF MS/MS modes. Leave set to 15.
Pre-Filter	2	The DC voltage applied to the pre and post filters of the quad. Set to 2V (DO NOT TUNE).
Ion Energy	0.2	The offset voltage applied to the quadrupole.

- Serves as ion guide in the single MS mode
- Served as mass filter in the tandem MS mode (precursor ion selection)

- Collision cell**
- Collision energy: 20-30 eV (collision energy is tuned by adjusting the voltage of focusing lenses, which accelerates ions into the collision cell)
  - Collision gas: argon
  - Pressure up to 10<sup>-4</sup> torr in CID
  - Collision RF values: RF offset is 150 for small molecules, and 400 for proteomics



**Quadrupole theory**

- 4 hyperbolic or cylindrical rods with the combination of dc (U) and rf voltages (V) applied on the rods
- Motion of ion in quadrupole fields:
 
$$a_x = 8zeU/mr^2\omega^2$$

$$q_x = 4zeV/mr^2\omega^2$$
- A constant a/q ratio will generate a scan line. Scanning this line (U, V varied) allows ions of differing m/z through the Q

**Quad profiles**

The Q will transmit over the region of ~4 times the set mass

**oa-TOF**

Ions are accelerated orthogonally (90°) by the pusher voltage

- Pusher: 1900 V
- Puller: 1400V
- Reflectron grid: 1700 V
- Reflectron: 1400V

- Data is acquired at speeds of up to 30,000 spectra per second
- Individual spectra are summed and a composite spectrum is recorded as raw data